



BEHRINGER

Medical gas devices



VACAVB area valve

Product description

The VAC AVB area valve is used to shut off VACUUM line pipes; it consists of a Ø 28 mm chemically nickel-plated copper pipe, a 1" ball valve and a G 1/4" F R threaded stub end for the installation of a pressure gauge or a pressure gauge with electrical contacts. Weld-on pipe inlet and outlet connections.

Normatives

UNI EN ISO 9001:2008 | UNI EN ISO 15001:2012

Components

- N. One Ø 28 mm chemically nickel-plated copper pipe.
- N. One 1" ball valve. N. 1 G1" 1/4 R - 1" R nipple.
- N. One G1/4" F R stub end for installing the pressure gauge.

Product codes

GAS	CODE	CATALOGUE NO.
VAC	HG096	1244969/R

Technical data

CODE	GAS	Ø PIPE
HG096	VAC	Ø28

Behringer S.r.l. via Gualco 4, 16165 Genova ITALY

T +39 010 8309103

F +39 010 8309183

e-mail info@behringer.it

C.F. e P.I. 01167840998