




BEHRINGER

Medical gas devices

Shut-off valve

Product description



Service shut-off valves are required to shut off each riser, branch and circuit of a medical gas and vacuum system. Each specific gas supply line must be shut off to complete maintenance or work downstream of the valve. Valve sizes range from 3/8" to 2" and the 3-piece weld-on pipe connectors at the outlets range from Ø 10 mm to Ø 54 mm. On the Ø 3/4" valves and larger, the 3-piece connectors are fitted with O-Rings to provide the best seal under pressure. All materials are oxygen compatible, supplied clean and free of any traces of oil, grease or dust.

Normatives

UNI EN ISO 7396-1:2013 | UNI EN ISO 9001:2008 | UNI EN ISO 15001:2012

Components

- N.1 ball valve with butterfly lever.
- N.Two 3-piece brass connectors

Product codes

VALVE SIZE	CODE
3/8"	H0670B
3/8"	H0671
1/2"	H0672
1/2"	H0672A
3/4"	H0674
1"	H0676
1" 1/4	H0677
1" 1/2	H0678
2"	H0679

Behringer S.r.l. via Gualco 4, 16165 Genova ITALY

T +39 010 8309103

F +39 010 8309183

e-mail info@behringer.it

C.F. e P.I. 01167840998

Technical data

CODE	VALVE SIZE	Ø PIPE	P. max
H0670B	3/8"	10 mm	63 bar
H0671	3/8"	12 mm	63 bar
H0672	1/2"	16 mm	50 bar
H0672A	1/2"	18 mm	50 bar
H0674	3/4"	22 mm	50 bar
H0676	1"	28 mm	45 bar
H0677	1" 1/4	35 mm	30 bar
H0678	1" 1/2	42 mm	25 bar
H0679	2"	54 mm	20 bar

Technical Images

DIAGRAMMA PRESSIONE/TEMPERATURA

