

Medical gas devices

UNI 9507 dual socket splitter adapter

Product description

The dual splitter adapter can be connected to one supply point so two UNI 9507 sockets can be used to dispense medical gas to two patients from one single supply point. All materials are oxygen compatible, supplied clean and free of any traces of oil, grease or dust.



Normatives

UNI EN ISO 9170-1:2008 | UNI EN ISO 9001:2008 | UNI 9507:2004 | UNI EN ISO 9170-2:2008 | UNI CEI EN ISO 13485:2004 | UNI EN ISO 15001:2012

Components

The dual splitter adapter body is made of CW614N chemically nickel-plated brass.

N.One UNI 9507 probe.

N.Two UNI 9507 sockets.

N.One PVC name ring.

Colored thermoplastic rubber ring.

Adhesive PVC gas label.

Weight 0.700 kg.

Product codes

GAS	CODE	CATALOGUE NO.
OXYGEN	HP059	158607/R
AIR 400	HP060	158607/R
VAC	HP061	158607/R
NITROUS OXIDE	HP062	158607/R