



# Medical gas devices

## AFNOR NF S 90-116 dual socket splitter adapter

### Product description

The dual splitter adapter can be connected to one supply point so two AFNOR NF S 90-116 sockets can be used to dispense medical gas to two patients from one single supply point. All materials are oxygen compatible, supplied clean and free of any traces of oil, grease or dust.



### Normatives

UNI EN ISO 9170-1:2008 | NF S 90-116 | UNI EN ISO 9001:2008 | UNI EN ISO 9170-2:2008 | UNI CEI EN ISO 13485:2004 | UNI EN ISO 15001:2012

### Components

The dual splitter adapter body is made of CW614N chemically nickel-plated brass.

- N. One AFNOR NF S 90-116 probe.
- N. Two AFNOR NF S 90-116 sockets.
- N. One PVC name ring.
- Colored thermoplastic rubber ring.
- Adhesive PVC gas label.
- Weight 0.450 kg.

### Product codes

GAS	CODE
OXYGEN	HP055
AIR 400	HP056
VAC	HP057
NITROUS OXIDE	HP058