



## UNI 9507 - AGSS EN ISO 9170-2 base block

### Product description

The UNI 9507 - AGSS EN ISO 9170-2 base block or base is made of CW614N chemically nickel-plated pressure die-cast brass with a Ø 16 mm welded copper inlet pipe.

### Normatives

UNI EN ISO 9001:2008 | UNI 9507:2004 | UNI EN ISO 9170-2:2008

### Components

Body made of CW614N chemically nickel-plated brass with a Ø 16 mm weld-on inlet pipe.

### Product codes

GAS	CODE	CATALOGUE NO.
AGSS	HP021A	159537/R

### Technical data

CODE	GAS
HP021A	AGSS

