



BEHRINGER

Medical gas devices



RVB collection jar

Product description

The collection jar is installed downstream of the RVB vacuum regulator to protect the suction system and regulator from impurities; it has a hose connector on the inlet and is connected to the vacuum regulator by a G 1/2" F R thread. The collection jar is made of polycarbonate and has a capacity of 150 cc or 300 cc. The top has an outlet with a hose connector.

Normatives

UNI EN ISO 9001:2008 | UNI CEI EN ISO 13485:2004 | UNI EN ISO 15001:2012

Components

Polycarbonate jar.

Yellow PVC top.

Moplen float.

One hose connector inlet.

N. One outlet with G 1/2" F R rotary nut. NBR O-ring.

Product codes

VOLUME	CONNECTION	CODE	CATALOGUE NO.
150 cc	G1/4" M R	HR081	217492/R
300 cc	G1/4" M R	HR082	217492/R

Technical data

CODE	VOLUME	CONNECTION
HR081	150 cc	G1/4" M R
HR082	300 cc	G1/4" M R

Behringer S.r.l. via Gualco 4, 16165 Genova ITALY

T +39 010 8309103

F +39 010 8309183

e-mail info@behringer.it

C.F. e P.I. 01167840998

