



BEHRINGER

Medical gas devices



AGSS ISO 9170-2

Product description

Anesthetic gas can be administered through the probe. The AGSS probe is ISO 9170-2 standard compliant.

Normatives

UNI EN ISO 7396-2:2007 | UNI CEI EN ISO 14971:2012 | UNI EN ISO 9170-2:2008 | UNI EN ISO 15001:2012

Components

Body made of CW614N chemically nickel-plated brass.
Specific label for each gas.

Product codes

CONNECTION	GAS	TYPE	CODE	CATALOGUE NO.
Hose connector Ø 13	AGSS	1 l	HI018	159659/R
Hose connector Ø 17	AGSS	1 l	HI018A	159659/R
Hose connector Ø 22	AGSS	1 l	HI019	159659/R
G3/8" F R	AGSS	1 l	HI019A	159659/R

Technical data

CODE	TYPE	CONNECTION
HI018	1 l	Hose connector Ø 13
HI018A	1 l	Hose connector Ø 17
HI019	1 l	Hose connector Ø 22
HI019A	1 l	G3/8" F R

Behringer S.r.l. via Gualco 4, 16165 Genova ITALY

T +39 010 8309103

F +39 010 8309183

e-mail info@behringer.it

C.F. e P.I. 01167840998

